

INSTRUCTION MANUAL

DNA Clean & Concentrator™-5

Catalog Nos. D4003T, D4003, D4004, D4013 & D4014

Highlights

- Quick, 2 minute recovery of ultra-pure DNA from PCR, enzymatic reactions, and other sources.
- DNA can be eluted in as little as 6 μl and is ideal for DNA ligation, sequencing, labeling, PCR, microarray, transfection, transformation, restriction digestion, etc.

Contents

Product Contents	1
Specifications	1
Product Description	2-3
Buffer Preparation	4
Protocol	4
Appendices ssDNA Purification Troubleshooting	
Ordering Information	7
List of Related Products	8-10

Page 1	
For Research Use Only	Ver. 1.2.1

Satisfaction of all Zymo Research products is guaranteed. If you should be dissatisfied with this product please call 1-888-882-9682.

Product Contents

DNA Clean & Concentrator ™-5 (Kit Size)	D4003T (10 Preps.)	D4003, D4013 (50 Preps.)	D4004, D4014 (200 Preps.)	Storage Temperature
DNA Binding Buffer	10 ml	50 ml	2 x 100 ml	Room Temp.
DNA Wash Buffer ¹	6 ml	6 ml	24 ml	Room Temp.
DNA Elution Buffer	1 ml	1 ml	4 ml	Room Temp.
Zymo-Spin™ Columns	10 Uncapped	50 D4003 – uncapped D4013 – capped	200 D4004 – uncapped D4014 – capped	Room Temp.
Collection Tubes	10	50	200	Room Temp.
Instruction Manual	1	1	1	-

Integrity of kit components is guaranteed for up to one year from date of purchase. Reagents are routinely tested on a lot-to-lot basis to ensure they provide the highest performance and reliability.

Specifications

- **DNA Purity** High-quality DNA ($A_{260}/A_{280} > 1.8$) ideal for ligation, sequencing, labeling, PCR, microarray, transfection, transformation, and restriction digestion procedures.
- DNA Size Limits From ~50 bp to 23 kb.
- DNA Recovery Typically, up to 5 μg total DNA per column can be eluted into as little as 6 μl of low salt DNA Elution Buffer or water. For DNA 50 bp to 10 kb, the recovery is 70-90%. For DNA 11 kb to 23 kb, the recovery is 50-70%.
- Sample Sources DNA from enzymatic reactions (e.g., PCR, restriction endonuclease digestions), plasmid preparations, and impure preparations.
- Product Detergent Tolerance ≤5% Triton X-100, ≤5% Tween-20, ≤5% Sarkosyl,
 ≤ 0.1% SDS.

Note: ™ Trademarks of Zymo Research Corporation. This product is for research use only and should only be used by trained professionals. It is not intended for use in diagnostic procedures. Some reagents included with this kit are irritants. Wear protective gloves and eye protection. Follow the safety guidelines and rules enacted by your research institution or facility. NanoDrop® is a registered trademark of NanoDrop Technologies, Inc.


¹ Ethanol must be added prior to use as indicated on the **DNA Wash Buffer** label. DNA Wash Buffer included with D4003S and D4003T is supplied ready-to-use and does not require the addition of ethanol prior to use.

Product Description


The <u>DNA Clean & Concentrator™-5</u> (DCC™-5) provides a hassle-free method for the rapid purification and concentration of high-quality DNA from PCR, endonuclease digestions, cell lysates, and other impure DNA preparations. It can also be used for post-RT cDNA clean-up and purification of sequencing-ready DNA from M13 phage. Simply add the specially formulated **DNA Binding Buffer** to your sample and transfer the mixture to the supplied **Zymo-Spin™ Column**. There is no need for organic denaturants or chloroform. Instead, the product features *Fast-Spin* column technology to yield DNA that is free of salts and contaminants in just 2 minutes. The purified DNA is ideal for DNA ligation, sequencing, labeling, PCR, microarray, transfection, transformation, and restriction digestion procedures.



Two minute **DCC-5™** procedure.



DNA samples, such as the PCR products shown here, can be efficiently purified and concentrated using the DNA Clean & Concentrator™-5.



Pure and Reliable Recovery with the DCCTM-5. Shown here is the recovery of 1 μg of 100 bp marker DNA eluted into 10 μl of water analyzed using a NanoDrop[®] spectrophotometer. The DNA Clean & Concentrator TM-5 consistently recovers > 90% of input DNA.

Available Formats

	DCC™-5	DCC™-25	DCC™-100	DCC™-500	Genomic DCC™	ZR-96 DCC™-5
Name	Zymo-Spin™ I & IC	Zymo-Spin™ II & IIC	Zymo-Spin™ V	Zymo-Spin™ VI	Zymo-Spin ™ IC-XL	Zymo-Spin™ I-96
Capacity	5 μg/ prep.	25 μg/ prep.	100 μg/ prep.	500 μg/ prep.	10 μg/ prep.	5 μg/ prep.
Elution Vol.	≥ 6 µl	≥ 25 µl	≥ 150 µl	≥ 2 ml	≥ 10 µl	≥ 10 µl
Cat. Nos.	D4003, D4013	D4005, D4033	D4029, D4030	D4031, D4032	D4010, D4011	D4023, D4024

Typical DCC™ Applications

Post-PCR DNA Clean-up	Efficient desalting of DNA with the removal of DNA polymerases, primers and free dNTPs.
DNA Clean-up From Enzymatic Reactions	Efficient desalting of DNA with the removal of modifying enzymes, RNA polymerases, ligases, kinases, nucleases, phosphatases, endonucleases, etc.
Post-Reverse Transcription (RT) & cDNA Clean-up	Efficiently purifies DNA following RT, either as a DNA/RNA complex or as single stranded cDNA following chemical hydrolysis of the RNA template.
Plasmid DNA Clean-up	Efficiently purifies plasmid DNA from "home-made" preparations of cell free lysates or from commercial kits. Plasmid DNA purified and concentrated using the DCC™ has proven an excellent substrate for high quality DNA sequencing.
Isotope and Dye Removal	Efficiently removes unincorporated fluorescent (<i>i.e.</i> , AMCA, FITC, BIO, DIG, Cy3, Cy5, FAM, <i>etc.</i>) and radiolabeled dNTP derivatives from DNA following <i>in vitro</i> labeling reactions.
Purification of M13 ssDNA	The DCC™ can be used for the rapid isolation of single stranded M13 phage DNA directly from phage-infected <i>E. coli</i> culture supernatant.

- ✓ For purification of short DNA or RNA oligonucleotides ≥ 16 nt, use the Oligo Clean & Concentrator (D4060, D4061).
- ✓ For ChIP (Chromatin Immunoprecipitation) sample cleanup, use the ChIP DNA Clean & Concentrator (D5201, D5205) for high quality DNA from any step in a standard ChIP protocol.
- ✓ For post-cycle sequencing samples, use the ZR Sequencing DNA Clean-up Kit (D4050, D4051) for dye blob elimination.
- ✓ For samples containing PCR inhibitors, use the *OneStep™ PCR Inhibitor Removal Kit* (D6030, D6035).

Selected Citations

Li, N. (2010). Whole genome DNA methylation analysis based on high throughput sequencing technology. *Methods*, *52* (3), 221-232. Lee, EJ. (2011). Targeted bisulfite sequencing by solution selection and massively parallel sequencing. *Nucleic Acids Research*, *39*(19), e127, doi:10.1093/nar/gkr598

Papageorgiou, EA. (2009). Sites of differential DNA methylation between placenta and peripheral blood. *Am J Pathol, 174* (5), 1609-1618. Ferguson, A.A. et al. (2009). Retrofitting ampicillin resistant vectors by recombination for use in generating *C. elegans* transgenic animals by bombardment. *Plasmid, 62,* 140-145.

Buffer Preparation

- ✓ <u>Before starting</u>: Add 24 ml 100% ethanol (26 ml 95% ethanol) to the 6 ml **DNA Wash Buffer** concentrate. Add 96 ml 100% ethanol (104 ml 95% ethanol) to the 24 ml **DNA Wash Buffer** concentrate.
- ✓ DNA Wash Buffer included with D4001S and D4001T is supplied ready-to-use and does not require the addition of ethanol prior to use.

For Assistance, please contact Zymo Research Technical Support at 1-888-882-9682 or e-mail tech@zymoresearch.com.

Protocol

All centrifugation steps should be performed between 10,000 - 16,000 x g.

1. In a 1.5 ml microcentrifuge tube, add 2-7 volumes of **DNA Binding Buffer** to each volume of DNA sample (see table below). Mix briefly by vortexing.

Application	DNA Binding Buffer: Sample	Example
Plasmid, genomic DNA (>2 kb)	2:1	200 µl : 100 µl
PCR product, DNA fragment	5 : 1	500 µl : 100 µl
ssDNA ¹ (e.g. cDNA, M13 phage)	7:1	700 µl : 100 µl

For efficient recovery of genomic or large DNA (> 20 kb to > 200 kb), use the **Genomic DNA Clean & Concentrator™ (Cat. Nos. D4010, D4011)**.

- 2. Transfer mixture to a provided **Zymo-Spin™ Column**² in a **Collection Tube**.
- 3. Centrifuge for 30 seconds. Discard the flow-through.
- 4. Add 200 μl **DNA Wash Buffer** to the column. Centrifuge for 30 seconds. Repeat the wash step.
- 5. Add ≥ 6 µl **DNA Elution Buffer**³ or water⁴ directly to the column matrix and incubate at room temperature for one minute. Transfer the column to a 1.5 ml microcentrifuge tube and centrifuge for 30 seconds to elute the DNA.

Ultra-pure DNA is now ready for use.

Notes:

- ¹ For ssDNA purification, see **Appendix A** on page 5.
- 2 The sample capacity of the column is 800 $\mu l.$ Therefore, it may be necessary to load and spin a column multiple times if a sample has a volume larger than 800 $\mu l.$
- ³ **DNA Elution Buffer**: 10 mM Tris-HCl, pH 8.5, 0.1 mM EDTA
- ⁴ Elution of DNA from the column is dependent on pH and temperature. If water is used, make sure the pH is >6.0. Waiting 1 minute prior to elution may improve the yield of larger (> 6 kb) DNA. For even larger DNA (> 10 kb), the total yield may be improved by eluting the DNA with 60-70 °C DNA Elution Buffer.

Appendix A: ssDNA Purification

cDNA clean-up

For clean-up of short cDNAs or ESTs (≥ 16 nt), we recommend the Oligo Clean & Concentrator ™ (Cat. Nos. D4060, D4061).

The DCC™ kit can be used to effectively clean and concentrate <u>cDNA</u> (> 500 nt) following reverse transcription (RT) in the presence/absence of fluorescent dyes. Unincorporated free nucleotides and fluorescent derivatives are efficiently removed using the DCC™, and the recovered cDNA may be used directly for microarray analysis, second-strand cDNA synthesis, or indirect labeling with a fluorescent dye such as NHS ester Cy3 or Cy5.

Hydrolysis

1. Add 10 μl 0.5 M EDTA and 10 μl 1 N NaOH to 50 μl of RT reaction.

The volumes of EDTA and NaOH should be scaled proportionally depending on the starting volume of the RT reaction.

2. Incubate at 65°C for 15 minutes.

Clean-up

 Add 490 µl (7 volumes) of **DNA Binding Buffer** to the hydrolysis reaction above. Mix well.

Neutralization (pH) following RNA hydrolysis is not necessary as the **DNA Binding Buffer** will effectively neutralize the NaOH added to the reaction.

2. Continue with Step 2 of the Protocol on page 4.

M13 phage ssDNA purification

- 1. Centrifuge phage-infected bacterial culture at 8,000 x q for 1 minute
- 2. Transfer 100 μl of phage-containing supernatant to a 1.5 ml microcentrifuge tube and add 700 μl (7 volumes) of **DNA Binding Buffer**. Mix briefly by vortexing.

Increased supernatant volumes may be processed by proportionally increasing the amount of **DNA Binding Buffer** added to the sample.

3. Continue with Step 2 of the Protocol on page 4.

Appendix B: Troubleshooting

Low Recovery

• Improperly Prepared/Stored DNA Wash Buffer

Make sure ethanol has been added to the **DNA Wash Buffer** concentrate. Cap the bottle tightly to prevent evaporation over time.

• Addition of DNA Elution Buffer

Add elution buffer directly to the column matrix and not to the walls of the column. Elution buffer requires contact with the matrix for at least 1 minute for large DNA \geq 10 kb.


• Incomplete Elution

- DNA elution is dependent on pH, temperature, and time. For large genomic DNA (≥ 50 kb), apply heated elution buffer (60-70 °C) and incubate for several minutes prior to elution.
- 2. Sequential elutions may be performed for quantitatively higher recovery but lower final DNA concentration. This is recommended for DNA ≥ 10 kb.

Low A₂₆₀/A₂₃₀ Ratios

• Column Tip Contaminated

When removing the column from the collection tube, be careful that the tip of the column does not come into contact with the flowthrough. Trace amounts of salt from the flowthrough can contaminate a sample resulting in low A_{260}/A_{230} ratios. Ethanol contamination from the flowthrough can also interfere with DNA elution. Zymo-SpinTM columns are designed for complete elution with no buffer retention or carryover (see below).



Following Clean-up with the DCC™, Multiple Bands Appear in an Agarose Gel

Acidification of DNA Loading Dye

Most loading dyes do not contain EDTA and will acidify (pH \leq 4) over time due to some microbial growth. This low pH is enough to cause DNA degradation. Therefore, if water is used to elute the DNA, 6X Loading Dye containing 1 mM EDTA is recommended.

Ordering Information

Product Description	Catalog No.	Kit Size (Preps.)
DNA Clean & Concentrator™-5 (for purification of up to 5 µg DNA per prep.) Supplied with uncapped columns	D4003T D4003 D4004	10 50 200
DNA Clean & Concentrator TM -5 (for purification of up to 5 μg DNA per prep.) Supplied with capped columns	D4013 D4014	50 200
ZR-96 DNA Clean & Concentrator [™] -5 (for 96-well purification of up to 5 µg DNA per well)	D4023 D4024	2 x 96 4 x 96
DNA Clean & Concentrator TM -25 (for purification of up to 25 μg DNA per prep.) Supplied with uncapped columns	D4005 D4006	50 200
DNA Clean & Concentrator TM -25 (for purification of up to 25 μg DNA per prep.) Supplied with capped columns	D4033 D4034	50 200
DNA Clean & Concentrator™-100 (for purification of up to 100 µg DNA per prep.)	D4029 D4030	25 50
DNA Clean & Concentrator TM -500 (for purification of up to 500 μg DNA per prep.)	D4031 D4032	10 20
Oligo Clean & Concentrator TM (for purification of up to 5 μg of oligonucleotides per prep.)	D4060 D4061	50 200
Genomic DNA Clean & Concentrator™ (for purification of up to 10 μg genomic DNA per prep.)	D4010 D4011	25 100

Refer to Page 3 for column design specifics in each kit.

For Individual Sale	Catalog No.	Size
DNA Binding Buffer	D4003-1-L D4004-1-L	50 ml 100 ml
DNA Wash Buffer (concentrate)	D4003-2-6 D4003-2-24	6 ml 24 ml
DNA Elution Buffer	D3004-4-1 D3004-4-4	1 ml 4 ml
Zymo-Spin™ I Columns (uncapped)	C1003-50 C1003-250	50 columns 250 columns
Zymo-Spin™ IC Columns (capped)	C1004-50 C1004-250	50 columns 250 columns
Collection Tubes	C1001-50 C1001-500 C1001-1000	50 tubes 500 tubes 1000 tubes

Popular Products From Zymo Research

Product	Description	Kit Size (Preps.)	Catalog No. (Format)			
	DNA Clean-up, Concentration & Recovery					
DNA Clean & Concentrator™-5	Clean and concentrate up to 5 µg DNA into ≥6 µl elution volume in as little as 2 minutes with no wash residue carryover.	50 200 50 200	D4003 (uncapped) D4004 (uncapped) D4013 (capped) D4014 (capped)			
DNA Clean & Concentrator™-25	Clean and concentrate 25 µg of DNA into ≥25 µl elution volume in as little as 2 minutes with no wash residue carryover.	50 200 50 200	D4005 (uncapped) D4006 (uncapped) D4033 (capped) D4034 (capped)			
ZR-96 DNA Clean & Concentrator™-5	Quick (30 minute), high throughput recovery of up to 5 μg pure DNA into 10-15 μl minimum elution volume allows for highly concentrated DNA.	2 x 96 4 x 96	D4023 D4024			
Genomic DNA Clean & Concentrator™	Quick (5 minute) clean-up of up to 10 µg high molecular weight DNA (≥ 20kb - 200 kb) from any enzymatic reaction or impure preparation without precipitations.	25 100	D4010 D4011			
Zymoclean™ Gel DNA Recovery Kit	Purify DNA from high and low-melting agarose gels in minutes.	50 200 50 200	D4001 (uncapped) D4002 (uncapped) D4007 (capped) D4008 (capped)			
ZR-96 Zymoclean™ Gel DNA Recovery Kit	High-throughput DNA purification from high and low-melting agarose gels.	2 x 96 4 x 96	D4021 D4022			
Zymoclean™ Large Fragment DNA Recovery Kit	Purify high molecular weight DNA (≥ 20 kb - 200 kb) from high and low-melting agarose gels in minutes.	25 100	D4045 D4046			
OneStep™ PCR Inhibitor Removal Kit	Fast, one step procedure for removal of PCR inhibitors such as polyphenolics, humic/fulvic acids, melanin, etc. for successful PCR and other downstream applications.	50 2 x 96	D6030 D6035			

Plasmid DNA Purification				
Zyppy™ Plasmid Miniprep Kit	Pellet-Free™ plasmid DNA purification in less than 10 minutes. Recover up to 25 μg DNA in as low as 30 μl.	50 100 400	D4036 D4019 D4020	
Zyppy™-96 Plasmid Miniprep	The fastest and simplest high-throughput method for plasmid purification. Magnetic bead format available for automated liquid handling platforms.	2 x 96 4 x 96 8 x 96 2 x 96 4 x 96 8 x 96	D4041 (spin plate) D4042 (spin plate) D4043 (spin plate) D4100 (magnetic bead) D4101 (magnetic bead) D4102 (magnetic bead)	
Zyppy™ Plasmid Midiprep Kit	Pellet-Free™ plasmid DNA purification in 15 minutes in a 150 µl minimum elution volume.	25 50	D4025 D4026	
ZR Plasmid MiniPrep™- <i>Classic</i>	Plasmid DNA purification in minutes: (alkaline lysis/spin column format for low 30 µl elution volume).	100 400 800	D4015 D4016 D4054	

Genomic DNA Purification				
<i>Quick-gDNA</i> ™ MiniPrep	Easy purification from whole blood, plasma, serum, body fluids, buffy coat, tissue, swabs or cultured cells ≥15 minutes without the use of Proteinase K or organic denaturants.	50/200 50/200	D3006/D3007 uncapped) D3024/D3025 (capped)	
ZR Genomic DNA™- Tissue MiniPrep	High quality DNA purification from solid tissues (e.g., tail snips, ear punches, adipose tissue, etc.), body fluids, cultured cells, buccal cells, FFPE tissues, hair, and other biological sources using Proteinase K and Fast.	50 200	D3050 D3051	
Environmental DNA Purification Kits	Unique BashingBead™ technology allows isolation of DNA from samples refractory to conventional lysis procedures including tough-to-lyse tissues, soil samples, feces, plants, seeds, insects, bacteria, yeast, filamentous fungi, unicellular and filamentous algae, and protozoa	Spin Column & 96-well Plate	Visit website for a comprehensive list	

RNA Purification				
RNA Clean & Concentrator™-5	Clean and concentrate up to 5 µg RNA into ≥6 µl elution volume in as little as 5 minutes with no wash residue carryover.	50 200	R1015 R1016	
Direct-Zol™ RNA MiniPrep	Quick, spin column purification of high-quality (DNA-free) total RNA <i>directly</i> from <i>TRI-Reagent</i> ® or similar acid-guanidinium-phenol based reagents (TRIzol®, RNAzol®, QIAzol®, TriPure, RNA-Bee <i>etc.</i>).	50 200	R2051 R2053	
ZR RNA MiniPrep	Rapid (15 minute) RNA isolation from a variety of sources using <i>Fast-Spin</i> column technology without the use of organic denaturants.	50 200	R1064 R1065	

Epigenetics Products From Zymo Research

Product	Description	Kit Size	Cat. No. (Format)
	Bisulfite Kits for DNA Methylation Detection		
EZ DNA Methylation™ Kit	For the conversion of unmethylated cytosines in DNA to uracil via the chemical-denaturation of DNA and a specially designed CT Conversion Reagent. Fast-Spin technology ensures ultra-pure, converted DNA for subsequent DNA methylation analysis. Magnetic bead format for automated liquid handling platforms.	50/200 Rxns. 2 x 96 Rxns. 2 x 96 Rxns. 4 x 96 Rxns.	D5001/D5002 (column) D5003 (shallow-well plate) D5004 (deep-well plate) D5040 (magnetic bead)
EZ DNA Methylation- Gold™ Kit	For the fast (3 hr.) conversion of unmethylated cytosines in DNA to uracil via heat/chemical-denaturation of DNA and a specially designed CT Conversion Reagent. Fast-Spin technology ensures ultra-pure, converted DNA for subsequent DNA methylation analysis. Magnetic bead format for automated liquid handling platforms.	50/200 Rxns. 2 x 96 Rxns. 2 x 96 Rxns. 4 x 96 Rxns.	D5005/D5006 (column) D5007 (shallow-well plate) D5008 (deep-well plate) D5042 (magnetic bead)
EZ DNA Methylation- Direct™ Kit	Simple and reliable DNA bisulfite conversion directly from blood, tissue (FFPE/LCM), and cells without the prerequisite for DNA purification in as little as 4-6 hrs. The increased sensitivity of this kit makes it possible to amplify bisulfite converted DNA from as few as 10 cells or 50 pg DNA. Magnetic bead format for adaptation to automated liquid handling platforms.	50/200 Rxns. 2 x 96 Rxns. 2 x 96 Rxns. 4 x 96 Rxns.	D5020/D5021 (column) D5022 (shallow-well plate) D5023 (deep-well plate) D5044 (magnetic bead)
EZ DNA Methylation- Lightning™ Kit	Complete bisulfite conversion in about an hour using a unique liquid format conversion reagent. Fast-Spin technology ensures ultra-pure, converted DNA for subsequent DNA methylation analysis. Magnetic bead format for automated liquid handling platforms.	50/200 Rxns. 2 x 96 Rxns. 2 x 96 Rxns. 4 x 96 Rxns.	D5030/D5031 (column) D5032 (shallow-well plate) D5033 (deep-well plate) D5046 (magnetic bead)
EZ DNA Methylation- Startup™ Kit	Designed for the first time user requiring a consolidated product to perform DNA methylation analysis. Includes technologies for sample processing, bisulfite treatment of DNA, and PCR amplification of "converted" DNA for methylation analysis.	1 Kit	D5024
	Methylated DNA Standards		1
Universal Methylated Human DNA Standard	Human (male) genomic DNA having all CpG sites methylated. To be used for the evaluation of bisulfite-mediated conversion of DNA. Supplied with a control primer set.	1 Set	D5011
Universal Methylated Mouse DNA Standard	Mouse (male) DNA having all CpG sites methylated. To be used for the evaluation of bisulfite-mediated conversion of DNA. Supplied with a control primer set.	1 Set	D5012
	Region-Specific DNA Methylation Screening		
<i>OneStep</i> qMethyl™ Kit	Single step real-time PCR procedure for bisulfite-free determination of DNA methylation status. Available without fluorescent dye for probe-based detection (Lite).	1 x 96 Rxns. 1 x 96 Rxns.	D5310 D5311 (Lite)
<i>OneStep</i> qMethyl™ Array	Premade 96-well assay for bisulfite-free determination of region-specific DNA methylation assessment in the promoter region of any one of the following prominent tumor suppressor genes: RASSF1, RARB, CDKN2A (p16), MGMT, or CCND2.	1 x 96 Rxns.	D5312
	Epigenetics Services		
Services for Methylated DNA	For more information, visit http://www.zymoresearch.com/services or inquire at services (A Analysis	@zymoresearch.co	<u>om</u> .
Simplify biomarker discovery	with our 5-mC Analysis platforms that combine Zymo's well-established bisulfite technologic tion analysis services available.	es with next-genera	ation sequencing for the most
Services for Hydroxymethyl Novel genome-wide 5-hmC ar sensitivity of 5-hmC detection	nalysis platform featuring cutting-edge 5-hmC DNA enrichment, library prep, and next-gene	ration sequencing	technologies to ensure the
	Hydroxymethylation Detection		
Quest 5-hmC™ DNA Enrichment Kit	Featuring J-base binding protein (JBP) for the specific enrichment of 5-hmC containing DNA, the consolidated workflow makes the procedure reliable for robust analysis of multiple samples.	25 Rxns. 50 Rxns.	D5420 D5421
Quest 5-hmC™ DNA ELISA Kit	Streamlined workflow for both the direct and relative quantitation of 5-hmC, in a global genomic context, with a robust colorimetric readout.	1 x 96 Rxns. 2 x 96 Rxns.	D5425 D5426
Anti-5- Hydroxymethylcytosine Polyclonal Antibody	Polyclonal antibody has been engineered to maximize sensitivity to low amounts of hydroxymethylated gDNA while minimizing crossreactivity with unmodified or methylated cytosine residues. The antibody is suitable for use in ELISA, IP, and immunohistochemical labeling.	50 μg 200 μg	A4001-50 A4001-200
DNA Degradase™ DNA Degradase Plus™	Whole genomic DNA can be treated with these enzyme cocktails for processing to individual nucleotides (Degradase™) or nucleosides (Degradase Plus™) for interrogation in chromatographic and spectroscopic methods including TLC, LC/MS, MALDI-TOF, and more.	500 U 2000 U 250 U 1000 U	E2016 E2017 E2020 E2021
	Other		
Zymo <i>Taq</i> ™ DNA Polymerase	Zymo Taq [™] "hot start" DNA Polymerase is specifically designed for the amplification of "difficult" DNA templates including: bisulfite-treated DNA for methylation detection. The product generates specific amplicons with little or no by-product formation.	50 Rxns. 200 Rxns	E2001/E2001 (system) E2003/E2004 (premix)
Methylated-DNA IP Kit	IP with a highly specific anti-5-methylcytosine monoclonal antibody. Designed for the enrichment of 5-methylcytosine-containing DNA from any pool of fragmented genomic DNA for use in genome-wide methylation analysis.	10 Rxns.	D5101




Premade Z-Competent[™] *E. coli* for 20 Second Transformations

(>108 transformants/µg DNA)

- √ NO Heat Shock!
- ✓ NO Lengthy Incubations!
- √ NO Outgrowth Procedures!
- ✓ NO Wait!!

Premade Z-Competent™ *E. coli* Cells

Product	Cat. No.	Size	
C600	T3015	10 x 100 μl aliquots (10 tubes)	
Zymo 5α	T3007	10 x 100 μl aliquots (10 tubes)	
(Same as DH5α)	T3009	96 x 50 μl aliquots (96-well plate)	
HB101	T3011	10 x 100 μl aliquots (10 tubes)	
ПВЮТ	T3013	96 x 50 μl aliquots (96-well plate)	
JM109	T3003	10 x 100 μl aliquots (10 tubes)	
JIVITOS	T3005	96 x 50 μl aliquots (96-well plate)	
TG1	T3017	10 x 100 μl aliquots (10 tubes)	
XJa Autolysis [™]	T3021	10 x 100 µl aliquots (10 tubes), 1 ml 500X L-Arabinose	
XJa(DE3) Autolysis™	T3031	10 x 100 µl aliquots (10 tubes), 1 ml 500X L-Arabinose	
XJb Autolysis [™]	T3041	10 x 100 µl aliquots (10 tubes), 1 ml 500X L-Arabinose	
XJb(DE3) Autolysis [™]	T3051	10 x 100 µl aliquots (10 tubes), 1 ml 500X L-Arabinose	



Make Your Own Z-Competent™ *E. coli* Cells

Product	Cat. No.	Size
Z-Competent [™] <i>E. coli</i> Transformation Kit (ZymoBroth [™] included)	T3001	up to 20 ml
Z-Competent [™] <i>E. coli</i> Transformation Buffer Set (ZymoBroth [™] not included)	T3002	up to 60 ml
Zuman Durath ™	M3015-100	100 ml
ZymoBroth [™]	M3015-500	500 ml